

QUALITY CONTROL CERTIFICATE

CICLOPIROX 50 µg Neo-Sensitabs

REF: 81112

LOT: 1601-1

Date of analysis: 2016-09

Date of expiry: 2021-08

Appearance: Light yellow, plane tablet, 9 mm in diameter, on both sides printed with a black identification code.
Complies.

Thickness: Mean of 10: 1.50-1.65 mm.
Complies.

Performance:

Test strains:

Candida albicans ATCC 64550


Candida albicans ATCC 10231

Product complies with specifications.

Uniformity: Inhibition zones of individual tablets vary less than ± 1 mm from mean.
Complies.

Microbiological quality: Complies.

Rosco Diagnostica A/S



Michela Venturi
Manager Quality Assurance