## **QUALITY CONTROL CERTIFICATE**

## CICLOPIROX 50 µg Neo-Sensitabs

REF: 81112N

LOT: 1401-1

Date of analysis: 07-2014

Date of expiry: 07-2019

Appearance: Light yellow, plane tablet, 9 mm in diameter,

on both sides printed with a black identification

code. Complies.

Thickness: Mean of 10: 1.50-1.65 mm.

Complies.

Performance:

Test strains:

Candida albicans ATCC 64550 Candida albicans ATCC 10231

Product complies with specifications.

<u>Uniformity:</u> Inhibition zones of individual tablets

vary less than  $\pm 1$  mm from mean.

Complies.

Microbiological quality: Complies.

ROSCO DIAGNOSTICA A/S

love B. Cares.

José B. Casals

Quality Assurance Manager