QUALITY CONTROL CERTIFICATE

CLOTRIMAZOLE 10 µg Neo-Sensitabs

REF: 81212N

LOT: 1202-1

Date of analysis: 09-2012

Date of expiry: 09-2017

Appearance: Light yellow, plane tablet, 9 mm in diameter,

on both sides printed with a black identification

code. Complies.

Thickness: Mean of 10: 1.50-1.65 mm.

Complies.

Performance:

Test strains:

Candida albicans ATCC 10231 Candida pseudotropicalis MH 2337

Product complies with specifications.

<u>Uniformity:</u> Inhibition zones of individual tablets

vary less than ± 1 mm from mean.

Complies.

Microbiological quality: Complies.

ROSCO DIAGNOSTICA A/S

José B. Casals

Quality Assurance Manager