

QUALITY CONTROL CERTIFICATE

CLOTRIMAZOLE 10 µg Neo-Sensitabs

REF: 81212

LOT: 1701-2

Date of analysis: 2017-09

Date of expiry: 2022-08

Appearance: Light yellow, plane tablet, 9 mm in diameter,
on both sides printed with a black identification
code.
Complies.

Thickness: Mean of 10: 1.50-1.65 mm.
Complies.

Performance:

Test strains:


Candida albicans ATCC 10231
Candida pseudotropicalis MIHE 2337

Product complies with specifications.

Uniformity: Inhibition zones of individual tablets
vary less than ± 1 mm from mean.
Complies.

Microbiological quality: Complies.

Rosco Diagnostica A/S



Michela Venturi
Manager Quality Assurance