QUALITY CONTROL CERTIFICATE

CLOTRIMAZOLE 10 µg Neo-Sensitabs

REF: 81212

LOT: 1701-3

Date of analysis: 2017-09

Date of expiry: 2022-08

<u>Appearance:</u> Light yellow, plane tablet, 9 mm in diameter,

on both sides printed with a black identification

code.
Complies.

Thickness: Mean of 10: 1.50-1.65 mm.

Complies.

Performance:

Test strains:

Candida albicans ATCC 10231 Candida pseudotropicalis MIHE 2337

Product complies with specifications.

<u>Uniformity:</u> Inhibition zones of individual tablets

vary less than ± 1 mm from mean.

Complies.

Microbiological quality: Complies.

Dewclalo Vento

Rosco Diagnostica A/S

Michela Venturi

Manager Quality Assurance