

## QUALITY CONTROL CERTIFICATE

### **ECONAZOLE 10 µg Neo-Sensitabs**

REF: 81312N

LOT: 1401-1

Date of analysis: 07-2014

Date of expiry: 07-2019

Appearance: Complies

Thickness: Mean of 10: 1.50-1.65 mm.  
Complies.

Performance:

Test strains:

*Candida albicans* ATCC 90028

*Candida parapsilosis* ATCC 22019

*Candida krusei* ATCC 6258

Product complies with specifications.

Uniformity: Inhibition zones of individual tablets  
vary less than  $\pm 1$  mm from mean.  
Complies.

Microbiological quality: Complies.

ROSCO DIAGNOSTICA A/S



---

Ole Skyggebjerg  
Quality Assurance Manager