QUALITY CONTROL CERTIFICATE

FLUOROCYTOSINE 1 µg Neo-Sensitabs

REF: 81512

LOT: 1601-2

Date of analysis: 2016-03

Date of expiry: 2021-02

Light yellow, plane tablet, 9 mm in diameter, Appearance:

on both sides printed with a black identification

code. Complies.

Mean of 10: 1.50-1.65 mm. Thickness:

Complies.

Performance:

Test strains:

Candida albicans ATCC 64548 Candida albicans ATCC 64550 Candida parapsilosis ATCC 22019

Candida krusei ATCC 6258

Product complies with specifications.

Uniformity: Inhibition zones of individual tablets

vary less than ± 1 mm from mean.

Complies.

Complies. Microbiological quality:

Dewclelo Vent

Rosco Diagnostica A/S

Michela Venturi

Manager Quality Assurance