## **QUALITY CONTROL CERTIFICATE**

## ITRACONAZOLE 10 µg Neo-Sensitabs

REF:	81812N
LOT:	1302-1
Date of analysis: Date of expiry:	04-2013 04-2018
<u>Appearance:</u>	Light yellow, plane tablet, 9 mm in diameter, on both sides printed with a black identification code. Complies.
Thickness:	Mean of 10: 1.50-1.65 mm. Complies.
Deufenneer	

Performance:

Test strains:

Candida parapsilosis ATCC 22019 Candida albicans ATCC 90028 Candida kruseii ATCC 6258

Product complies with specifications.

Uniformity:

Inhibition zones of individual tablets vary less than  $\pm 1$  mm from mean. Complies.

Microbiological quality: Complies.

ROSCO DIAGNOSTICA A/S

ve B. Cars.

José B. Casals Quality Assurance Manager