## QUALITY CONTROL CERTIFICATE

## ITRACONAZOLE 10 µg Neo-Sensitabs

REF:	81812N
LOT:	1401-2
Date of analysis: Date of expiry:	2014-08 2019-07
<u>Appearance:</u>	Light yellow, plane tablet, 9 mm in diameter, on both sides printed with a black identification code. Complies.
<u>Thickness:</u>	Mean of 10: 1.50-1.65 mm. Complies.
Dorformanca	

Performance:

Test strains:

Candida parapsilosis ATCC 22019 Candida albicans ATCC 90028 Candida kruseii ATCC 6258

Product complies with specifications.

Uniformity:

Inhibition zones of individual tablets vary less than  $\pm 1$  mm from mean. Complies.

Microbiological quality: Complies.

ROSCO DIAGNOSTICA A/S

ve B. Cars.

José B. Casals Quality Assurance Manager