## **QUALITY CONTROL CERTIFICATE**

## ITRACONAZOLE 10 μg Neo-Sensitabs

REF: 81812N

LOT: 1401-1

Date of analysis: 2014-08 Date of expiry: 2019-07

<u>Appearance:</u> Light yellow, plane tablet, 9 mm in diameter,

on both sides printed with a black identification

code.
Complies.

Thickness: Mean of 10: 1.50-1.65 mm.

Complies.

Performance:

Test strains:

Candida parapsilosis ATCC 22019 Candida albicans ATCC 90028 Candida kruseii ATCC 6258

Product complies with specifications.

Uniformity: Inhibition zones of individual tablets

vary less than  $\pm 1$  mm from mean.

Complies.

Microbiological quality: Complies.

ROSCO DIAGNOSTICA A/S

José B. Casals

Quality Assurance Manager