QUALITY CONTROL CERTIFICATE

ITRACONAZOLE 10 µg Neo-Sensitabs

REF: 81812

LOT: 2001-1

Date of analysis: 2020-05

Date of expiry: 2025-03

<u>Appearance:</u> Light yellow, plane tablet, 9 mm in diameter,

on both sides printed with a black identification

code. Complies.

Thickness: Mean of 10: 1.50-1.65 mm.

Complies.

Performance:

Test strains:

Candida parapsilosis ATCC 22019 Candida albicans ATCC 90028 Candida kruseii ATCC 6258

Product complies with specifications.

Uniformity: Inhibition zones of individual tablets

vary less than ± 1 mm from mean.

Complies.

Microbiological quality: Complies.

Rosco Diagnostica A/S

Peter Samuelsen

Manager Quality Assurance