## QUALITY CONTROL CERTIFICATE

## NYSTATIN 50 µg Neo-Sensitabs

REF:	82212N
LOT:	1401-1
Date of analysis:	06-2014
Date of expiry:	06-2019
<u>Appearance:</u>	Light yellow, plane tablet, 9 mm in diameter, on both sides printed with a black identification code. Complies.
Thickness:	Mean of 10: 1.50-1.65 mm. Complies.
Performance:	

Test strains:

*Candida albicans* ATCC 64548 *Candida albicans* ATCC 10231

Product complies with specifications.

Uniformity:Inhibition zones of individual tablets<br/>vary less than  $\pm 1$  mm from mean.<br/>Complies.

Microbiological quality: Complies.

ROSCO DIAGNOSTICA A/S

hose B. Cars.

José B. Casals Quality Assurance Manager