QUALITY CONTROL CERTIFICATE

NYSTATIN 50 µg Neo-Sensitabs

REF: 82212

LOT: 1701-3

Date of analysis: 2017-04

Date of expiry: 2022-04

<u>Appearance:</u> Light yellow, plane tablet, 9 mm in diameter,

on both sides printed with a black identification

code. Complies.

Thickness: Mean of 10: 1.50-1.65 mm.

Complies.

Performance:

Test strains:

Candida albicans ATCC 64548 Candida albicans ATCC 10231

Product complies with specifications.

Uniformity: Inhibition zones of individual tablets

vary less than ± 1 mm from mean.

Complies.

Microbiological quality: Complies.

Dewclelo Vent

Rosco Diagnostica A/S

Michela Venturi

Manager Quality Assurance