QUALITY CONTROL CERTIFICATE

VORICONAZOLE 1 µg Neo-Sensitabs

REF: 82312N

LOT: 1101-1

Date of analysis: 03-2011

Date of expiry: 03-2016

Appearance: Light yellow, plane tablet, 9 mm in diameter,

on both sides printed with a black identification

code. Complies.

Thickness: Mean of 10: 1.50-1.65 mm.

Complies.

Performance:

Test strains:

Candida albicans ATCC 90028 Candida parapsilesis ATCC 22019 Candida krusei ATCC 6258

Product complies with specifications.

<u>Uniformity:</u> Inhibition zones of individual tablets

vary less than ± 1 mm from mean.

Complies.

Loss on drying: 0.46 %

Complies.

Microbiological quality: Complies.

ROSCO DIAGNOSTICA A/S

Mikkel Pandrup Duer-Jensen Manager Quality Assurance