QUALITY CONTROL CERTIFICATE

CASPOFUNGIN 5 µg Neo-Sensitabs

REF: 82412

LOT: 2001-1

Date of analysis: 2020-04

Date of expiry: 2022-10

Appearance: Light yellow, plane tablet, 9 mm in diameter,

on both sides printed with a black identification

code.
Complies.

Thickness: Mean of 10: 1.50-1.65 mm.

Complies.

Performance:

Test strains:

Candida albicans ATCC 90028 Candida parapsilosis ATCC 22019 Candida krusei ATCC 6258

Product complies with specifications.

<u>Uniformity:</u> Inhibition zones of individual tablets

vary less than ± 1 mm from mean.

Complies.

•

Microbiological quality: Complies.

Rosco Diagnostica A/S

Peter Samuelsen

Manager Quality Assurance