## **QUALITY CONTROL CERTIFICATE**

## FLUCONAZOLE 25 µg Neo-Sensitabs

REF: 82512N

LOT: 1302-1

Date of analysis: 05-2013

Date of expiry: 05-2018

<u>Appearance:</u> Light yellow, plane tablet, 9 mm in diameter,

on both sides printed with a black identification

code. Complies.

Thickness: Mean of 10: 1.50-1.65 mm.

Complies.

Performance:

Test strains:

Candida Paraposilosis ATCC 22019 Candida albicans ATCC 90028

Product complies with specifications.

<u>Uniformity:</u> Inhibition zones of individual tablets

vary less than  $\pm 1$  mm from mean.

Complies.

.

Microbiological quality: Complies.

ROSCO DIAGNOSTICA A/S

José B. Casals

Quality Assurance Manager