QUALITY CONTROL CERTIFICATE

FLUCONAZOLE 25 µg Neo-Sensitabs

REF: 82512N

LOT: 1501-1

Date of analysis: 2015-07

Date of expiry: 2020-07

Appearance: Light yellow, plane tablet, 9 mm in diameter,

on both sides printed with a black identification

code. Complies.

Thickness: Mean of 10: 1.50-1.65 mm.

Complies.

Performance:

Test strains:

Candida Paraposilosis ATCC 22019 Candida albicans ATCC 90028

Product complies with specifications.

<u>Uniformity:</u> Inhibition zones of individual tablets

vary less than ± 1 mm from mean.

Complies.

.

Microbiological quality: Complies.

Rosco Diagnostica a/s

Dorrit Reumert

Quality Assurance Manager