

## Certificate of Analysis

**Bruker Bacterial Test Standard** Product Name:

Part No.

8255343

LOT No.:

0000320889

Use by<sup>1)</sup>:

2021-11-30

**Date of Production:** 

2019-12-10

QC release date:

2020-01-20

**Storage Conditions:** Store at -18 °C or below immediately on arrival.

<u>Test</u>	<b>Specification</b>	LOT Results
Verification of the calibration points by MALDI-TOF MS	Within a mass range of ± 900ppm	Passed
Relative intensity of reference proteins No. 3 by MALDI-TOF MS	$0.20 \pm 0.05$	Passed
Relative intensity of reference proteins No. 4 by MALDI-TOF MS	$0.22 \pm 0.05$	Passed
MALDI Biotyper identification	Classification Result Report lists "Escherichia coli" as best match.	Passed
Remarks	$\log(\text{score})$ Value $\geq 2.2$ of the best match	Passed

## Remarks:

1) BTS is shipped at ambient temperature but must be stored at -18°C or below immediately on arrival. The expiry date on the package is valid for the enclosed BTS when stored at -18°C or below on arrival. Frozen, dissolved Bruker BTS can be stored for up to 5 months at -18°C or below. See instruction for use at www.bruker.com/care,

Bruker Daltonik GmbH, Germany Fahrenheitstr. 4 – 28359 Bremen phone +49-421-2205-0

Dr. Markus Kostrzewa, SAFETY OFFICER