

• Certificate of Analysis

Product Name : Bruker Bacterial Test Standard

Part No. 8255343

LOT No.: 7030220013

Use by¹⁾: 2022-08-31

Date of Production: 2020-09-16

QC release date: 2020-09-25

Storage Conditions : Store at -18 °C or below immediately on arrival.

| <u>Test</u> | <u>Specification</u> | <u>LOT Results</u> |
|--|--|--------------------|
| Verification of the calibration points by MALDI-TOF MS | Within a mass range of ± 900 ppm | Passed |
| Relative intensity of reference proteins No. 3 by MALDI-TOF MS | 0.20 ± 0.05 | Passed |
| Relative intensity of reference proteins No. 4 by MALDI-TOF MS | 0.22 ± 0.05 | Passed |
| MALDI Biotyper identification | Classification Result Report lists "Escherichia coli" as best match. | Passed |
| | log(score) Value ≥ 2.2 of the best match | Passed |

Remarks:

1) BTS is shipped at ambient temperature but must be stored at -18°C or below immediately on arrival. The expiry date on the package is valid for the enclosed BTS when stored at -18°C or below on arrival. Frozen, dissolved Bruker BTS can be stored for up to 5 months at -18°C or below. See instruction for use at www.bruker.com/care.

Bruker Daltonik GmbH, Germany
Fahrenheitstr. 4 – 28359 Bremen
phone +49-421-2205-0



Dr. Markus Kostrzewa, SAFETY OFFICER