

Certificate of Analysis

Product Name: Bruker Bacterial Test Standard

 Part No.
 8255343

 LOT No.:
 7030220017

 Use by¹¹:
 2022-10-31

 Date of Production:
 2020-11-26

 QC release date:
 2021-01-05

Storage Conditions: Store at -18 °C or below immediately on arrival.

<u>Test</u>	Specification	LOT Results
Verification of the calibration points by MALDI-TOF MS	Within a mass range of ± 900ppm	Passed
Relative intensity of reference proteins No. 3 by MALDI-TOF MS	0.20 ± 0.05	Passed
Relative intensity of reference proteins No. 4 by MALDI-TOF MS	0.22 ± 0.05	Passed
MALDI Biotyper identification	Classification Result Report lists "Escherichia coli" as best match.	Passed
	$log(score)$ Value ≥ 2.2 of the best match	Passed

Remarks:

1) BTS is shipped at ambient temperature but must be stored at -18°C or below immediately on arrival. The expiry date on the package is valid for the enclosed BTS when stored at -18°C or below on arrival. Frozen, dissolved Bruker BTS can be stored for up to 5 months at -18°C or below. See instruction for use at www.bruker.com/care.

Bruker Daltonik GmbH, Germany Fahrenheitstr. 4 – 28359 Bremen phone +49-421-2205-0

Bruker Daltonik GmbH Fahrenheitstraße A IS-28359 Bremen

Dr. Markus Kostrzewa, SAFETY OFFICER