

Certificate of Analysis

Product Name:

Bruker Matrix HCCA, portioned

Part No .:

8255344

LOT No .:

0000251625

Use by¹⁾:

2018-07

Date of Production:

2016-08-12

Storage Conditions:

Store at 2°C to 8°C on arrival.

Test:

Specification

LOT Results

1. Amount of HCCA content. Sampling in accordance to DIN ISO 3951-1.

AQL-value=1.0 Inspection level II U=2.8 mg // L=2.2 mg Passed

2. MALDI Biotyper identification of a 1:20 diluted BTS.

Sampling in accordance to DIN ISO

Classification Result Report lists "Escherichia coli" as best match. Passed

2859-1.

AQL-value=1.0

Passed

Inspection level S4

 $L=2.000 \log(\text{score})$ of the best match

3. Blank value analysis of HCCAportioned by MALDI Biotyper. Sampling in accordance to DIN ISO

Classification Result Report lists "no reliable identification".

Passed

2859-1.

AQL-value=1.0 Inspection level S4 $U=1.699 \log(score)$

Passed

Remarks:

1) Bruker HCCA is shipped at ambient temperature. The expiry date on the package is valid for the product when stored in a refrigerator at 2°C to 8°C on arrival.

Bruker Daltonik GmbH, Germany Fahrenheitstr. 4 – 28359 Bremen phone +49-421-2205-0

QC release date: 01-09-2016

Hannah Tegtmeyer - QA Departmen

Not for use in diagnostic procedures. For research use only.