

Certificate of Analysis

Product Name:

Bruker Matrix HCCA, portioned

Part No.

8255344

LOT No.:

0000242534

Use by¹⁾:

2018-02

Date of Production:

2016-03-07

Storage Conditions:

Store at 2°C to 8°C on arrival.

<u>Test</u>	Specification	LOT Results
Differential weighing	Tolerance: ± 0.3 mg	Passed
Functional Check against Protein Calibration Standard I (Bruker)	Signal/Noise: ≥ 100	Passed
	Resolution: ≥ 250	Passed

Remarks:

1) Bruker HCCA is shipped at ambient temperature. The expiry date on the package is valid for the product when stored in a refrigerator at 2°C to 8°C on arrival. Dissolved Bruker HCCA is stable at controlled room temperature (20–25°C) for up to one week.

Bruker Daltonik GmbH, Germany Fahrenheitstr. 4 – 28359 Bremen phone +49-421-2205-0

QC release date: 08/03/2016

Certified Quality Quality Assurance

Not for use in diagnostic procedures. For research use only.