

LOT

## Certificate of Analysis

Product Name:

Bruker Matrix HCCA, portioned

Part No.:

8255344

LOT No.:

0000278004

Use by<sup>1)</sup>:

2019-08

Date of Production:

2017-09-12

Storage Conditions:

Store at 2°C to 8°C on arrival.

Test:	<b>Specification</b>	Results
1. Amount of HCCA content.	AQL-value=1.0	Passed
Sampling in accordance to DIN ISO	Inspection level II	
3951-1.	U=2.8 mg // L=2.2 mg	
2. MALDI Biotyper identification	Classification Result Report lists	Passed
of a 1:20 diluted BTS.	"Escherichia coli" as best match.	
Sampling in accordance to DIN ISO		
2859-1.	AQL-value=1.0	Passed
	Inspection level S4	
	L= 2.000 log(score) of the best match	
3. Blank value analysis of HCCA-	Classification Result Report lists "no	Passed
portioned by MALDI Biotyper.	reliable identification".	
Sampling in accordance to DIN ISO		
2859-1.	AQL-value=1.0	Passed
	Inspection level S4	

## Remarks:

1) Bruker HCCA is shipped at ambient temperature. The expiry date on the package is valid for the product when stored in a refrigerator at 2°C to 8°C on arrival.

U= 1.699 log(score)

Bruker Daltonik GmbH, Germany Fahrenheitstr. 4 – 28359 Bremen phone +49-421-2205-0

QC release date: 27-09-2017

Fahrenheidstraße 4

D-28359 Bremen
Tel: +49 (0) 421 - 22 05 - 0

r.M. Kostrzewa, SAFETY OFFICER

Not for use in diagnostic procedures. For research use only.