

Certificate of Analysis

Product Name:

Bruker Matrix HCCA, portioned

Part No .:

8255344

LOT No .:

0000282924

Use by1):

2020-01

Date of Production:

2018-02-14

Storage Conditions:

Store at 2°C to 8°C on arrival.

Test:

Specification

LOT Results

Passed

1. Amount of HCCA content. Sampling in accordance to DIN ISO 3951-1.

AQL-value=1.0 Inspection level II U=2.8 mg // L=2.2 mg

AOL-value=1.0

Passed

2. MALDI Biotyper identification of a 1:20 diluted BTS.

Sampling in accordance to DIN ISO

Classification Result Report lists "Escherichia coli" as best match.

2859-1.

Inspection level S4 L= 2.000 log(score) of the best match Passed

3. Blank value analysis of HCCAportioned by MALDI Biotyper. Sampling in accordance to DIN ISO 2859-1.

Classification Result Report lists "no reliable identification".

Passed

AQL-value=1.0 Inspection level S4 U= 1.699 log(score)

Passed

Remarks:

1) Bruker HCCA is shipped at ambient temperature. The expiry date on the package is valid for the product when stored in a refrigerator at 2°C to 8°C on arrival.

Bruker Daltonik GmbH, Germany Fahrenheitstr. 4 – 28359 Bremen phone +49-421-2205-0

ahrenheitstraße 4 -28359 Bremen Tel: +49 (0) 421 - 22 05 - 0

Bruker Daltonik GmbH

Dr.M. Kostrzewa, SAFETY OFFICER

QC release date: 21-02-2018

Not for use in diagnostic procedures. For research use only.