

Certificate of Analysis

Product Name:

Bruker Matrix HCCA, portioned

Part No .:

8255344

LOT No .:

0000293411

Use by1):

2020-08

Date of Production:

2018-09-11

Storage Conditions:

Store at 2°C to 8°C on arrival.

Test:

Specification

LOT Results

1. Amount of HCCA content. Sampling in accordance to DIN ISO AQL-value=1.0 Inspection level II U=2.8 mg // L=2.2 mg Passed

3951-1. 2. MALDI Biotyper identification

of a 1:20 diluted BTS. Sampling in accordance to DIN ISO

2859-1.

Classification Result Report lists "Escherichia coli" as best match.

Passed

AQL-value=1.0

Passed

Inspection level S4

L= 2.000 log(score) of the best match

3. Blank value analysis of HCCAportioned by MALDI Biotyper. Sampling in accordance to DIN ISO

Classification Result Report lists "no reliable identification".

Passed

2859-1.

AQL-value=1.0 Inspection level S4 $U=1.699 \log(score)$

Passed

Remarks:

1) Bruker HCCA is shipped at ambient temperature. The expiry date on the package is valid for the product when stored in a refrigerator at 2°C to 8°C on arrival.

Bruker Daltonik GmbH, Germany Fahrenheitstr. 4 – 28359 Bremen phone +49-421-2205-0

Tel: +49 (0) 421 - 22 05 - 0

Dr.M. Kostrzewa, SAFETY OFFICER

Bruker Daltonik GmbH

QC release date: 08-10-2018

Not for use in diagnostic procedures. For research use only.