

Certificate of Analysis

Product Name: Bruker Matrix HCCA, portioned

Part No.: 8255344

LOT No.: 0000300602

Use by1): 2020-10

Date of Production: 2018-11-14

Storage Conditions: Store at 2°C to 8°C on arrival.

LOT Test: **Specification** Results 1. Amount of HCCA content. AQL-value=1.0 Passed Sampling in accordance to DIN ISO Inspection level II 3951-1. U=2.8 mg // L=2.2 mgClassification Result Report lists 2. MALDI Biotyper identification Passed "Escherichia coli" as best match. of a 1:20 diluted BTS. Sampling in accordance to DIN ISO 2859-1. AOL-value=1.0 Passed Inspection level S4 L= 2.000 log(score) of the best match 3. Blank value analysis of HCCA-Classification Result Report lists "no Passed reliable identification". portioned by MALDI Biotyper. Sampling in accordance to DIN ISO AOL-value=1.0 2859-1. Passed

Remarks:

1) Bruker HCCA is shipped at ambient temperature. The expiry date on the package is valid for the product when stored in a refrigerator at 2°C to 8°C on arrival.

Inspection level S4 U= 1.699 log(score)

Bruker Daltonik GmbH, Germany Fahrenheitstr. 4 – 28359 Bremen phone +49-421-2205-0

QC release date: 30-11-2018

ruker Daltonik GmbH Arenheitstraße/4 Tel: +49 (0) 421 - 22 05 - 0 Dr.M. Kostrzewa, SAFETY OFFICER

Not for use in diagnostic procedures. For research use only.