

Passed

## Certificate of Analysis

Product Name: Bruker Matrix HCCA, portioned

Part No.: 8255344

LOT No .: 7020120079

Use by1): 2022-08

Date of Production: 2020-09-16

Storage Conditions: Store at 2°C to 8°C on arrival.

**LOT** Test: Specification Results

1. Amount of HCCA content. AQL-value=1.0 Sampling in accordance to DIN ISO Inspection level II 3951-1. U=2.8 mg // L=2.2 mg

Passed 2. MALDI Biotyper identification Classification Result Report lists "Escherichia coli" as best match. of a 1:20 diluted BTS.

Sampling in accordance to DIN ISO 2859-1. AQL-value=1.0 Passed

Inspection level S4  $L=2.000 \log(\text{score})$  of the best match

3. Blank value analysis of HCCA-Passed portioned by MALDI Biotyper. reliable identification". Sampling in accordance to DIN ISO 2859-1. AQL-value=1.0 Passed

Classification Result Report lists "no

Inspection level S4  $U=1.699 \log(\text{score})$ 

Remarks:

1) Bruker HCCA is shipped at ambient temperature. The expiry date on the package is valid for the product when stored in a refrigerator at 2°C to 8°C on arrival.

Bruker Daltonik GmbH, Germany Fahrenheitstr. 4 – 28359 Bremen phone +49-421-2205-0

QC release date: 25-09-2020

Bruker Daltonik GmbH Fahrenheitstraße 4 D-28359 Bremen Tel: +49 (0) 421 AFETY OFFICER

Not for use in diagnostic procedures. For research use only.