

Certificate of Analysis

Product Name: Bruker Matrix HCCA, portioned

Part No.:

8255344

LOT No .:

6020121006

Use by1):

2023-06

Date of Production:

2021-07-01

Storage Conditions:

Store at 2°C to 8°C on arrival.

LOT Test: Specification Results

1. Amount of HCCA content. Sampling in accordance to DIN ISO 3951-1.

AQL-value=1.0 Inspection level II U=2.8 mg // L=2.2 mg Passed

2. MALDI Biotyper identification of a 1:20 diluted BTS.

Sampling in accordance to DIN ISO 2859-1.

Classification Result Report lists "Escherichia coli" as best match. Passed

AOL-value=1.0

Passed

Inspection level S4

L= 2.000 log(score) of the best match

3. Blank value analysis of HCCAportioned by MALDI Biotyper. Sampling in accordance to DIN ISO 2859-1.

Classification Result Report lists "no reliable identification".

Passed

AQL-value=1.0 Inspection level S4 $U=1.699 \log(\text{score})$

Passed

Remarks:

1) Bruker HCCA is shipped at ambient temperature. The expiry date on the package is valid for the product when stored in a refrigerator at 2°C to 8°C on arrival.

Bruker Daltonics GmbH & Co. KG, Germany Fahrenheitstr. 4 – 28359 Bremen phone +49-421-2205-0

Bruker Daltonics GmbH & Co. KG Fahrenheitstr. 4 D-28359 Bremen Phone + 49 (0)421-2205-0 ex +49 (0)423 AFETYOFFICER

QC release date: 14-07-2021

Not for use in diagnostic procedures. For research use only.