

Passed

Passed

Certificate of Analysis

Product Name: Bruker Matrix HCCA, portioned

Part No.: 8255344

LOT No.: 6020121048

Use by¹⁾: 2023-10

Date of Production: 2021-11-08

Storage Conditions: Store at 2°C to 8°C on arrival.

Test: Specification LOT Results

1. Amount of HCCA content.AQL-value=1.0Sampling in accordance to DIN ISOInspection level II3951-1.U=2.8 mg // L=2.2 mg

2. MALDI Biotyper identification Classification Result Report lists "Escherichia coli" as best match.

Sampling in accordance to DIN ISO
2859-1. AQL-value=1.0 Passed
Inspection level S4

L= 2.000 log(score) of the best match

3. Blank value analysis of HCCAportioned by MALDI Biotyper.
Sampling in accordance to DIN ISO

2859-1.

AQL-value=1.0
Inspection level \$4.

Inspection level S4 U= 1.699 log(score)

Remarks:

1) Bruker HCCA is shipped at ambient temperature. The expiry date on the package is valid for the product when stored in a refrigerator at 2°C to 8°C on arrival.

Bruker Daltonics GmbH & Co. KG, Germany Fahrenheitstr. 4 – 28359 Bremen phone +49-421-2205-0

QC release date: 17-11-2021

Not for use in diagnostic procedures. For research use only.

Phone + 49 (0)421-2205-0 Dr.M. Kostrzewa, SAFETOY OFFICER

Fahrenheitst.4 D-28359 Bremen

Bruker Daltonics GmbH & Co. KG C