QUALITY CONTROL CERTIFICATE

Posaconazole 5 µg Neo-Sensitabs

REF: 82612N

LOT: 1301-1

Date of analysis: 10-2013

Date of expiry: 09-2018

<u>Appearance:</u> Light yellow, plane tablet, 9 mm in diameter,

on both sides printed with a black identification

code. Complies.

Thickness: Mean of 10: 1.50-1.65 mm.

Complies.

Performance:

Test strains:

C Albicans ATCC 90028 C Krusei ATCC 6258

C Parapsilosis ATCC 22019

Product complies with specifications.

<u>Uniformity:</u> Inhibition zones of individual tablets

vary less than ± 1 mm from mean.

Complies.

•

Microbiological quality: Complies.

ROSCO DIAGNOSTICA A/S

José B. Casals

Quality Assurance Manager