QUALITY CONTROL CERTIFICATE

Posaconazole 5 µg Neo-Sensitabs

REF: 82612

LOT: 1601-2

Date of analysis: 2016-04

Date of expiry: 2021-04

Appearance: Light yellow, plane tablet, 9 mm in diameter,

on both sides printed with a black identification

code. Complies.

Thickness: Mean of 10: 1.50-1.65 mm.

Complies.

Performance:

Test strains:

C Albicans ATCC 90028 C Krusei ATCC 6258

C Parapsilosis ATCC 22019

Product complies with specifications.

<u>Uniformity:</u> Inhibition zones of individual tablets

vary less than ± 1 mm from mean.

Complies.

.

Microbiological quality: Complies.

Vendela Venta

Rosco Diagnostica A/S

Michela Venturi

Manager Quality Assurance