

QUALITY CONTROL CERTIFICATE

Posaconazole 5 µg Neo-Sensitabs

REF: 82612

LOT: 210009

Date of analysis: 2021-03

Date of expiry: 2026-03

Appearance: Light yellow, plane tablet, 9 mm in diameter,
on both sides printed with a black identification
code.
Complies.

Thickness: Mean of 10: 1.50-1.65 mm.
Complies.

Performance:

Test strains:

C Albicans ATCC 90028
C Krusei ATCC 6258
C Parapsilosis ATCC 22019

Product complies with specifications.

Uniformity: Inhibition zones of individual tablets
vary less than ± 1 mm from mean.
Complies.

Microbiological quality: Complies.

Rosco Diagnostica A/S



Peter Samuelson
Manager Quality Assurance