

## Certificate of Analysis

Product Name:

**IVD Bacterial Test Standard** 

Part No.

8290190

LOT No.:

0000320837

Use by<sup>1)</sup>:

2020-05-31

**Date of Production:** 

2019-06-19

QC release date:

2019-07-17

**Storage Conditions:** 

Store at -18 °C or below immediately on arrival.

<u>Test</u>	<b>Specification</b>	LOT Results
Verification of the calibration points by MALDI-TOF MS	Within a mass range of ± 900ppm	Passed
Relative intensity of reference proteins No. 3 by MALDI-TOF MS	$0.20 \pm 0.05$	Passed
Relative intensity of reference proteins No. 4 by MALDI-TOF MS	$0.22 \pm 0.05$	Passed
MALDI Biotyper identification	Classification Result Report lists "Escherichia coli" as best match.	Passed
Remarks	log(score) Value ≥ 2.2 of the best match	Passed

## Remarks

1) IVD BTS is shipped at ambient temperature but must be stored at -18°C or below immediately on arrival. The expiry date on the package is valid for the enclosed IVD BTS when stored at -18°C or below on arrival. Frozen, dissolved Bruker IVD BTS can be stored for up to 5 months at -18°C or below. See instruction for use at <a href="https://www.bruker.com/care">www.bruker.com/care</a>.

Bruker Daltonik GmbH, Germany Fahrenheitstr. 4 – 28359 Bremen phone +49-421-2205-0

Bruker Daltonik GmbH

5ahrenheitstraße 4

D 28359 Gremen

Tel: +49 (0) 421 - 22 05 - 0

Tel: +49(6) 421 - 22 05 - 0 Dr. Markus Kostrzewa, SAFETY OFFICER