

Certificate of Analysis

Product Name:

IVD Bacterial Test Standard

Part No.

8290190

LOT No.:

7030321001

Use by¹⁾:

2021-12-31

Date of Production:

2021-01-04

QC release date:

2021-01-13

Storage Conditions:

Store at -18 °C or below immediately on arrival.

<u>Test</u>	Specification	LOT Results
Verification of the calibration points by MALDI-TOF MS	Within a mass range of ± 900ppm	Passed
Relative intensity of reference proteins No. 3 by MALDI-TOF MS	0.20 ± 0.05	Passed
Relative intensity of reference proteins No. 4 by MALDI-TOF MS	0.22 ± 0.05	Passed
MALDI Biotyper identification	Classification Result Report lists "Escherichia coli" as best match.	Passed
	log(score) Value ≥ 2.2 of the best	Passed

Remarks:

1) IVD BTS is shipped at ambient temperature but must be stored at -18°C or below immediately on arrival. The expiry date on the package is valid for the enclosed IVD BTS when stored at -18°C or below on arrival. Frozen, dissolved Bruker IVD BTS can be stored for up to 5 months at -18°C or below. See instruction for use at www.bruker.com/care.

Bruker Daltonik GmbH, Germany Fahrenheitstr. 4 – 28359 Bremen phone +49-421-2205-0

Bruker Daltovik GmbH

Dr/Markus Kostrzewa, SAEE TY OFFICER