

Quality control certificate



Nutrient Broth N°2

N° Cat: 1251

N° Lote: 909022

Re-test Date: 9/2023

Formula in g/l

Casein peptone	10	Beef extract	10
Sodium chloride	5		

Quality control

Solubility	Appearance	Color of the dehydrated medium	Color of the prepared medium	Final pH (25°)	Lote pH
w/o rests	Fine powder	Beige	Amber, slightly opalescent	7,5 ± 0,2	7,5

Microbiological test

Incubation conditions: (35±2 °C / 18-48 h)

Microrganisms	Specification
Staphylococcus epidermidis ATCC 12228	Good growth
Streptococcus pyogenes ATCC 12344	Good growth
Enterobacter aerogenes ATCC 13048	Good growth
Escherichia coli ATCC 25922	Good growth
Salmonella typhi ATCC 6539	Good growth

This certificate is an electronic copy of the certificate available in our laboratory and does not require signature.