Quality control certificate



VL Anaerobic Blood Agar

N° Cat: 1485 N° Lote: 104214 Re-test Date: 4/2024

Formula in g/l

<u> </u>			
Glucose	2	Bacteriological agar	17,6
L-Cysteine hydrochloride	0,4	Peptone	10
Sodium chloride	5	Yeast extract	5
Concentrated Ox Extract	2		· · · · · · · · · · · · · · · · · · ·

Quality control

Solubility	Appearance	Color of the dehydrated medium	Color of the prepared medium	Final pH (25°)	Lote pH
w/o rests	Fine powder	Light toasted	Light amber. With blood: red cherry	6.9 ± 0.2	6.8

Microbiological test

Incubation conditions: (35±2 °C / 24-48 h) under anaerobic conditions

Microorganisms	Specification
Clostridium sporogenes ATCC 3485	Good growth

This certificate is an electronic copy of the certificate available in our laboratory and does not require signature.