

Quality control certificate



Bordetella Supplement

N° Cat: 6015

N° Lote: 320151

Re-test Date: 1/2023



Formula per vial

Cephalexin (mg)	20
-----------------	----

Quality control

Solubility	Appearance	Color of the dehydrated medium	Color of the prepared medium	Final pH (25°)	Lote pH
w/o rests	Lyophilized tablet	N/A	Transparent	N/A	0,0

Microbiological test

Incubation conditions: (35±2 °C / 48-72 h).

Microorganisms	Specification	Characteristic reaction
Bordetella parapertussis ATCC 15311	Good growth	Gamma hemolysis
Bordetella bronchiseptica ATCC 4617	Good growth	Gamma hemolysis
Bordetella pertussis ATCC 8467	Good growth	Beta hemolysis

This certificate is an electronic copy of the certificate available in our laboratory and does not require signature.