

SAFE® CAP 8 CHANNEL HH

Semi-Automated 8-Channel Capper/Decapper

- Semi-automated hand-held capper/decapper with ergonomic design
- Reliably decaps and caps eight tubes in parallel
- Stable stand for storage of device and contamination-free parking of screw caps
- Eight individual motors for precise application of 7 Ncm torque for tube closure
- Power-saving mode for longer operating time in wireless mode

Item Number

SAFE® 96 XT:

SAFE® 96 IT:

CDC-8CH-XT-2

CDC-8CH-IT-2

Order Information



SAFE® CAP 8 CHANNEL HH