

SAFE[®] CAP MULTI CAP DD

Automated Single Column/Row Capper/Decapper for SBS 96/48/24

- Fully-automated premium capper/decapper
- Easily swapped decapping heads to decap tubes in 24-, 48- and 96-format tube racks on a single device
- Process a full rack of tubes or row of tubes in portrait or landscape format within one device
- Decap only tubes needed — all rows, selected rows or columns, or partial racks
- Can be operated as an integrated and standalone device simultaneously
- Eliminates risk of cross contamination by not moving over opened tubes
- Minimize the time a tube is open using the optional Row Loop Mode — only one row is processed at a time by holding caps after decapping and immediately recapping
- Built-in Secure Mode ensures an optimal seal during capping to eliminate cross threading
- Touchscreen panel provides easy navigation through the system
- Turntable allows rack to be turned after decapping for better positioning in pipetting workflows
- Dimensions: 600 mm x 380 mm x 440 mm

Item Number

SAFE[®] 96 XT:

SAFE[®] 96 IT:

SAFE[®] 48 XT:

SAFE[®] 24 XT:

CDC-MC-12CH-XT-DD

CDC-MC-12CH-IT-DD

CDC-MC-8CH-XT-DD

CDC-MC-6CH-XT-DD

Order Information



SAFE[®] CAP MULTI CAP DD