



www.pro-lab.com

CERTIFICATE OF ANALYSIS

Product Name: IMMERSION OIL	Lot Number: 88598
Product Code: PL.396	Expiration Date: 2019 05

Additional Information

Description: Viscous liquid.	Concentrate Lot number: n/a Production date: 2014/7/11	Storage conditions: 15-25°C Precautions: refer to MSDS
------------------------------	---	---

Quality Control Testing

QC Results:
 Appearance: Pale straw coloured viscous liquid
 Refractive index at 20°C: 1.518
 Weight per ml: 0.93g
 Viscosity at 20°C: 1000-1500 cSt
 Colour < 50 Hazen Units
 ABBE Number 44

Test results:
 n/a

Comments:
 We hereby certify that samples of the above Lot have been tested in accordance with current specifications and that the above is a true record of the results obtained.

Release

14/07/11	<i>mawr .</i>
Date	Quality Control

Canada ☎ 800 268 2341 Fax 905 731 0206
 20 Mural Street, Unit #4, Richmond Hill, Ontario,
 L4B 1K3
 support@pro-lab.com

U.S.A. ☎ 800 522 7740 Fax 800 332 0450
 21 Cypress Blvd., Suite 1070, Round Rock,
 Texas 78665-1034
 ussupport@pro-lab.com

U.K. ☎ 0151 353 1613 Fax 0151 353 1614
 3 Bassendale Road, Bromborough, Wirral,
 Merseyside CH62 3QL
 uksupport@pro-lab.com