



www.pro-lab.com

CERTIFICATE OF ANALYSIS

Product Name: SAFRANIN	Lot Number: K30849
Product Code: PL.7012	Expiration Date: 2021 11 30

Additional Information

Description: Deep red solution with slight odour.	Concentrate Lot number: K30790 Production date: 2020/12/09	Storage conditions: 15-25°C Precautions: refer to SDS
---	---	--

Quality Control Testing

Challenged with test organisms:
Escherichia coli NCTC 10418
Oxford *Staphylococcus aureus* NCTC 6571
Haemolytic *Streptococcus* Group A NCTC 8198

Test results:
Conforms to performance specifications when used in standard Grams stain method, giving crisp staining of Gram positive and Gram negative bacteria with a clean counterstained background.

Comments:
We hereby certify that samples of the above Lot have been tested in accordance with current specifications and that the above is a true record of the results obtained.

Release

2020/12/10	
Date	Quality Control

Canada ☎ 800 268 2341 Fax 905 731 0206
20 Mural Street, Unit #4, Richmond Hill, Ontario,
L4B 1K3
support@pro-lab.com

U.S.A. ☎ 800 522 7740 Fax 800 332 0450
21 Cypress Blvd., Suite 1070, Round Rock,
Texas 78665-1034
ussupport@pro-lab.com

U.K. ☎ 0151 353 1613 Fax 0151 353 1614
3 Bassendale Road, Bromborough, Wirral,
Merseyside CH62 3QL
uksupport@pro-lab.com