



www.pro-lab.com

### CERTIFICATE OF ANALYSIS

<b>Product Name:</b> MALACHITE GREEN	<b>Lot Number:</b> K32994
<b>Product Code:</b> PL.8007	<b>Expiration Date:</b> 2023 01 31

#### Additional Information

Description: Dark green solution.	Concentrate Lot number: K32945 Production date: 2022/01/18	Storage conditions: 15-25°C Precautions: refer to SDS
-----------------------------------	---	--

#### Quality Control Testing

Challenged with test organisms:  
*Mycobacterium scrofulaceum* NCTC 10803  
*Escherichia coli* NCTC 12241

Test results:  
 Conforms to performance specifications when used in standard Ziehl Neelson stain method, giving crisp red staining of Acid/AlcoholFast bacilli with a clean counterstained background.

Comments:  
 We hereby certify that samples of the above Lot have been tested in accordance with current specifications and that the above is a true record of the results obtained.

#### Release

2022/01/24	
Date	Quality Control

Canada ☎ 800 268 2341 Fax 905 731 0206  
 20 Mural Street, Unit #4, Richmond Hill, Ontario,  
 L4B 1K3  
 support@pro-lab.com

U.S.A. ☎ (512) 832 9145 Fax (512) 832 6424  
 21 Cypress Blvd., Suite 1155, Round Rock,  
 Texas 78665-1046  
 support@pro-lab.us

U.K. ☎ 0151 353 1613 Fax 0151 353 1614  
 3 Bassendale Road, Bromborough, Wirral,  
 Merseyside CH62 3QL  
 uksupport@pro-lab.com