

CERTIFICATE OF ANALYSIS

Product Name: RABBIT COAGULASE PLASMA	Lot Number: 180446
Product Code: PL.850-3, PL.850-10, PL.850-30	Expiration Date: 2021 04

Additional Information		
Description: Rabbit plasma containing disodium EDTA as anticoagulant.	Pack size: 30ml Production date: 2018/04/23	Storage conditions: 2-8°C Precautions: refer to SDS

Quality Control Testing


Raw Materials and Components:
All incoming chemicals, raw materials and components are quality controlled against predetermined specifications to ensure reproducibility.

Production:
Production is carried out against a Product Manufacturing Procedure incorporating Standard Operating Procedures and Good Manufacturing Practices.

Finished Product:
Lot to lot comparisons are carried out to ensure that the appearance of the material, as well as the specificity and sensitivity, remains constant.

Quality Assurance
Testing against panels of *S.aureus* and *S.epidermidis* are carried out against predetermined parameters. The above lot passed all tests.

Stability and Shelf Life of Lyophilised Product
Recommended storage temperature:2-8°C. Storage at lower temperatures has no detrimental effect. Extended stability trials at 4°C, 25°C and 37°C have shown that the product has a minimum shelf life of 2 years at 4°C and at least one month at 25°C. Reconstituted material, using sterile distilled water, should not be stored longer than necessary and should never exceed the recommendations of the product literature.

Release	
2018/04/23	
Date	Quality Control

Canada ☎ 800 268 2341 Fax 905 731 0206
20 Mural Street, Unit #4, Richmond Hill, Ontario,
L4B 1K3
support@pro-lab.com

U.S.A. ☎ 800 522 7740 Fax 800 332 0450
21 Cypress Blvd., Suite 1070, Round Rock,
Texas 78665-1034
ussupport@pro-lab.com

U.K. ☎ 0151 353 1613 Fax 0151 353 1614
3 Bassendale Road, Bromborough, Wirral,
Merseyside CH62 3QL
uksupport@pro-lab.com