



www.pro-lab.com

CERTIFICATE OF ANALYSIS

Product Name: RABBIT COAGULASE PLASMA	Lot Number: 210437B
Product Code: PL.850-10	Expiration Date: 2024 04 30

Additional Information

Description: Rabbit plasma containing disodium EDTA as anticoagulant.	Pack size: 10ml Production date: 2021/05/10	Storage conditions: 2-8°C Precautions: refer to SDS
---	--	--

Quality Control Testing

Raw Materials and Components:

All incoming chemicals, raw materials and components are quality controlled against predetermined specifications to ensure reproducibility.

Production:

Production is carried out against a Product Manufacturing Procedure incorporating Standard Operating Procedures and Good Manufacturing Practices.

Finished Product:

Lot to lot comparisons are carried out to ensure that the appearance of the material, as well as the specificity and sensitivity, remains constant.


Quality Assurance

Testing against S.aureus NCTC 12981 / ATCC 25923, S.aureus NCTC 6571 / ATCC 9144 and S.epidermidis NCTC 13360 / ATCC 12228 is carried out against predetermined parameters. The above lot passed all tests.

Stability and Shelf Life of Lyophilised Product

Recommended storage temperature:2-8°C. Storage at lower temperatures has no detrimental effect. Extended stability trials at 4°C, 25°C and 37°C have shown that the product has a minimum shelf life of 2 years at 4°C and at least one month at 25°C. Reconstituted material, using sterile distilled water, should not be stored longer than necessary and should never exceed the recommendations of the product literature.

Release

2021/05/10	
Date	Quality Control

Canada ☎ 800 268 2341 Fax 905 731 0206
20 Mural Street, Unit #4, Richmond Hill, Ontario,
L4B 1K3
support@pro-lab.com

U.S.A. ☎ (512) 832 9145 Fax (512) 832 6424
21 Cypress Blvd., Suite 1155, Round Rock,
Texas 78665-1046
support@pro-lab.us

U.K. ☎ 0151 353 1613 Fax 0151 353 1614
3 Bassendale Road, Bromborough, Wirral,
Merseyside CH62 3QL
uksupport@pro-lab.com