

QUALITY CONTROL CERTIFICATE

OPTOCHIN 10 µg Diagnostic Tablets

REF: 44211

LOT: 1901-2

Date of analysis: 2019-04

Date of expiry: 2024-04

Appearance: Plane tablet, 6 mm in diameter.
Complies.

Average weight: 110 mg.
Complies.

Performance:

Test strains:

Streptococcus pneumoniae ATCC 49619

Streptococcus bovis ATCC 15351

Streptococcus mitis ATCC 17693

Product complies with specification.

Rosco Diagnostica A/S



Peter Samuelson
Manager Quality Assurance