

QUALITY CONTROL CERTIFICATE

BETA-GLUCURONIDASE Diagnostic Tablets PGUA

REF: 50611
LOT: 1103-1
Date of analysis: 12-2011
Date of expiry: 12-2016
Appearance: Plane tablet, 6 mm in diameter.
Complies.
Average weight: 90 mg.
Complies.

Performance:

Test strains:

E. coli ATCC 25922
Klebsiella pneumoniae ATCC 13883
Staphylococcus simulans ATCC 27848
S. haemolyticus ATCC 29970

Product complies with specification.

Loss on drying: 0.58 %

ROSCO DIAGNOSTICA A/S



Ole Skyggebjerg
Quality Assurance Manager