

## QUALITY CONTROL CERTIFICATE

### **MALTOSE Diagnostic Tablets**

REF: 52911  
LOT: 1201-1  
Date of analysis: 09-2012  
Date of expiry: 08-2016  
Appearance: Plane tablet, 6 mm in diameter.  
Complies.  
Average weight: 90 mg.  
Complies.


### Performance:

Test strains:

*E. coli* ATCC 25922  
*Morganella morganii* ATCC 25830  
*Staphylococcus aureus* ATCC 25923  
*Staphylococcus capitis* 3477b

Product complies with specification.

ROSCO DIAGNOSTICA A/S



José B. Casals  
Quality Assurance Manager