

## QUALITY CONTROL CERTIFICATE

### **Total Metallo-beta-lactamasas Confirm Kit**

REF: 98016

LOT: 1902-1

Date of expiry: 2021-12

Appearance: Light yellow, plane tablets, 9 mm in diameter,  
on both sides printed with a black identification  
code.  
Complies.

Thickness: Mean of 10: 1.50-1.65 mm.  
Complies.

Performance:

Test strains	Results
<i>E. coli</i> ATCC 25922	Negative
<i>K. pneumoniae</i> ATCC BAA - 2146	MBL Positive
<i>K. pneumoniae</i> NCTC 13438	MBL Positive

Product complies with specifications.

Uniformity: Inhibition zones of individual tablets  
vary less than  $\pm 1$  mm from mean.  
Complies.

Microbiological quality: Complies.

Rosco Diagnostica A/S



Peter Samuelson  
Manager Quality Assurance