

QUALITY CONTROL CERTIFICATE

Detection of COLISTIN resistance, Neo-Sensitabs

REF: 98018

LOT: 1301-1

Date of analysis: 05-2013

Date of expiry: 05-2018

Appearance: Light yellow, plane tablets, 9 mm in diameter,
on both sides printed with a black identification
code.
Complies.

Thickness: Mean of 10: 1.50-1.65 mm.
Complies.

Quality Control:
Pseudomonas aeruginosa ATCC 27853: 20-29 mm
E. coli ATCC 27853 : 30-38 mm

Product complies with specifications.

Uniformity: Inhibition zones of individual tablets
vary less than ± 1 mm from mean.
Complies.

Microbiological quality: Complies.

ROSCO DIAGNOSTICA A/S



José B. Casals
Quality Assurance Manager