

QUALITY CONTROL CERTIFICATE

KPC/MBL in *P. aeruginosa*/Acinetobacter Confirm Kit, Neo-Sensitabs

REF: 98020
LOT: 1301-1
Date of analysis: 12-2012, 11-2012, 06-2013, 12-2012,06-2013
Date of expiry: 2015-11

Appearance: Light yellow, plane tablets, 9 mm in diameter, on both sides printed with a black identification code.
Complies.

Thickness: Mean of 10: 1.50-1.65 mm.
Complies.

Performance:

Test strains:

Test strains	Results
<i>P. aeruginosa</i> ATTC 10145	MBL Positive
<i>P. aeruginosa</i> CCUG 59626	MBL Positive
<i>K. pneumoniae</i> CCUG 56233	KPC Positive
<i>K. pneumoniae</i> NCTC 13438	KPC Positive

Product complies with specifications.

Uniformity: Inhibition zones of individual tablets vary less than ± 1 mm from mean.
Complies.

Microbiological quality: Complies.

ROSCO DIAGNOSTICA A/S



José B. Casals
Quality Assurance Manager